DUKWOO MEDICAL

PRODUCT CATALOG

About us

- About DUKWOO
- History
- Quality

Products


- IV Catheter
- Intervention
- Centesis and Drainage
- Women's Healthcare
- General Purpose Catheter














Quality

For the purpose of quality management, we not only were awarded ISO, but also try to meet quality system requirements exceeding those of CE, FDA and KGMP.







Certificate of GMP Certificate of ISO Certificate of CE



Type: Straight & Blister package Sterilization: EO gas sterilization Non-toxic and non-pyrogenic

| MODEL | SPECIFICATION (Guage) | CATHETER LENGTH (inch) (mm) | | COLOR CODE |
|------------------|--------------------------|-----------------------------|-----|------------|
| 3S-Cath+ 14G 51 | | 2 | 51 | |
| 3S-Cath+ 14G 80 | | 3 | 80 | 0044165 |
| 3S-Cath+ 14G 100 | - 14 | 4 | 100 | ORANGE |
| 3S-Cath+ 14G 150 | | 6 | 150 | |
| 3S-Cath+ 16G 45 | | 1 3/4 | 45 | |
| 3S-Cath+ 16G 51 | | 2 | 51 | |
| 3S-Cath+ 16G 80 | 16 | 3 | 80 | GREY |
| 3S-Cath+ 16G 100 | | 4 | 100 | |
| 3S-Cath+ 16G 150 | | 6 | 150 | |
| 3S-Cath+ 18G 32 | | 1 1/4 | 32 | |
| 3S-Cath+ 18G 45 | | 1 3/4 | 45 | |
| 3S-Cath+ 18G 51 | 10 | 2 | 51 | CDEEN |
| 3S-Cath+ 18G 80 | - 18 | 3 | 80 | GREEN |
| 3S-Cath+ 18G 100 | | 4 | 100 | |
| 3S-Cath+ 18G 150 | | 6 | 150 | |
| 3S-Cath+ 20G 25 | | 1 | 25 | |
| 3S-Cath+ 20G 32 | 20 | 1 1/4 | 32 | RED |
| 3S-Cath+ 20G 51 | | 2 | 51 | |
| 3S-Cath+ 22G 19 | 22 | 3/4 | 19 | BLUE |
| 3S-Cath+ 22G 25 |] | 1 | 25 | DLUE |
| 3S-Cath+ 24G 19 | 24 | 3/4 | 19 | YELLOW |

Sheath-Pia Introducer

for interventional access



Sheath Introducer from Dukwoo Medical is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries.

Sheath

- Radiopaque material sheath / dilator enhances visualization.
- Longer tapered sheath and dilator improves vascular access.
- Reduce the vessel tissue trauma.

Guidewire

- Soft and flexible tip allows to access easily into the vessel.
- Reduce the vessel tissue trauma.

Components:

Sheath Introducer 1. Dilator 1. Guide Wire 1. Access Needle 1. Syringe 1. These components may be packaged together or separately.

| | MODEL SPECIFICATION | | | GUIDI | WIRE | | I.V.CA | THETER | | COLOR |
|---------------------|---------------------|-------------|--------------|----------------|---------|----------|-----------|----------------|---------|--------|
| MODEL | MODEL (Fr) | LENGTH (mm) | OD (inch) | LENGTH (mm) | TYPE | TIP | OD (G) | LENGTH (mm) | SYRINGE | CODE |
| Sheath-Pia-A410 | 4 | 100 | X | X | X | X | X | X | X | RED |
| Sheath-Pia-A510 | 5 | 100 | Х | X | X | Х | X | Х | X | GREY |
| Sheath-Pia-A610 | 6 | 100 | X | X | X | X | X | Х | Х | GREEN |
| Sheath-Pia-A710 | 7 | 100 | X | X | X | X | X | X | X | ORANGE |
| Sheath-Pia-A810 | 8 | 100 | X | X | X | X | X | X | X | BLUE |
| Sheath-Pia-E410-35S | 4 | 100 | 0.035 | 450 | Plastic | Straight | 18 | 51 | 0 | RED |
| Sheath-Pia-E510-35S | 5 | 100 | 0.035 | 450 | Plastic | Straight | 18 | 51 | 0 | GREY |
| Sheath-Pia-E610-35S | 6 | 100 | 0.035 | 450 | Plastic | Straight | 18 | 51 | 0 | GREEN |
| Sheath-Pia-E710-35S | 7 | 100 | 0.035 | 450 | Plastic | Straight | 18 | 51 | 0 | ORANGE |
| Sheath-Pia-E810-35S | 8 | 100 | 0.035 | 450 | Plastic | Straight | 18 | 51 | 0 | BLUE |
| Sheath-Pia-H57 | 5 | 70 | X | X | X | X | X | X | Х | GREY |
| Sheath-Pia-H67 | 6 | 70 | X | X | X | × | X | X | X | GREEN |
| Sheath-Pia-H77 | 7 | 70 | X | X | X | X | X | Х | X | ORANGE |

3S-Micro Introducer



Micro Introducer from Dukwoo Medical is intended to provide access and facilitate the percutaneous introduction of various devices.

- Sheath / Dilator
- Stiffer dilator increases the durability of the device.
- Longer tapered sheath and dilator improves vascular access.
- Radio opaque material sheath and dilator enhances visualization.

Components:

Micro Introducer Sheath 1. (inner dilator with outer sheath) 0.018" x 450mm Guide Wire 1.

Access Needle 1.

These components may be packaged together or separately.

Guidewire

- Nitinol with Gold distal tip.
- Radiopaque distal tip enhances visualization.
- Soft and flexible tip allows to access easily into the vessel.
- Nitinol wire reduces kinking into the vessel.

| MODEL | SPECIFICATION | CATHETER | NEEDLE | NEEDLE | GUIDEWIRE | | |
|-------------------------|---------------|----------------|--------------|----------------|--------------|----------------|----------|
| MODEL | (Fr) | LENGTH (mm) | GUAGE (G) | LENGTH (mm) | OD (inch) | LENGTH (mm) | TYPE |
| Sheath-Pia-M410 | | | Х | Х | X | Х | Х |
| Sheath-Pia-M410-18T-ang | 4 | 100 | 18 | 70 | 0.018 | 450 | angled |
| Sheath-Pia-M410-18T-S | | | 10 | 70 | 0.018 | 430 | Straight |
| Sheath-Pia-M510 | | | X | X | × | × | X |
| Sheath-Pia-M510-18T-ang | 5 | 100 | 18 | 70 | 0.018 | 450 | angled |
| Sheath-Pia-M510-18T-S | | | 18 | 70 | 70 0.018 | 450 - | Straight |

3S-Percutaneous Access Introducer 5Fr



3S-Percutaneous Access Introducer from Dukwoo Medical is percutaneous introduction of biliary catheters, a 5F sheath assembly with stiffening stylet.

This introducer allows for the introduction of 0.018" guide wire in percutaneous procedures.

Components:

Introducer Sheath 1. (inner stylet with outer sheath) 0.018" x 700mm Guide Wire 1.

Chiba Needle 1.

These components may be packaged together or separately.

Available in 5 French (5F x 300mm)

| MODEL | SPECIFICATION | SPECIFICATION CATHETER | | NEEDLE | GUIDEWIRE | | | |
|----------------------|---------------|------------------------|--------------|--------|--------------|-------------|--------|--|
| INIODEL | (Fr) | LENGTH (mm) | GUAGE (G) | | OD (inch) | LENGTH (mm) | TYPE | |
| N530-S | | | X | X | X | X | X | |
| N530-22G-18T-S-70ang | 5 | 270 | 22 | 150 | 0.018 | 700 | angled | |
| N530-21G-18T-S-70ang | | | 21 | 150 | 0.018 | 700 | angled | |

3S-Percutaneous Access Introducer 6Fr



3S-Percutaneous Access Introducer from Dukwoo Medical is percutaneous introduction of biliary catheters, a 6F sheath assembly with stiffening stylet. This introducer allows for the introduction of 0.018" guide wire in percutaneous procedures. Radiopaque marker band on the sheath

- Sheath / Dilator / Stiffening cannula
- Stiff Cannula increases the durability of the device.
- Longer tapered sheath ease introduction and placement.
- Radiopague material sheath and dilator enhances visualization.
- Very lubricious and low friction during insertion.
- Minimize the vessel tissue trauma and patient pain.
- Components:

Introducer Sheath 1. (stiffening cannula, dilator and sheath) 0.018" x 600mm Guide Wire 1.

Chiba Needle 1.

These components may be packaged together or separately.

- Guidewire
- Nitinol with Gold distal tip.
- Radiopaque distal tip enhances visualization.
- Soft and flexible tip allows to access easily into the vessel.
- Nitinol wire reduces kinking into the vessel.

| MODEL | SPECIFICATION | SPECIFICATION CATHETER | | NEEDLE | | | |
|--------------------|---------------|------------------------|-----------------------|--------|--------------|-------------|--------|
| MODEL | (Fr) | LENGTH (mm) | GUAGE LENGTH (G) (mm) | | OD (inch) | LENGTH (mm) | TYPE |
| N620-S-M | | | Х | X | X | × | X |
| N620-22G-18T-60ang | 6 | 185 | 22 | 150 | 0.018 | 700 | angled |
| N620-21G-18T-60ang | | | 21 | 150 | 0.018 | 700 | angled |



3S-A.D.Cath from Dukwoo Medical is intended to drainage ascites.

Components:

Ascites Drainage Catheter 1. (with hemostasis valve and side arm tube) Extension Tube 1. (1,060mm)

These components may be packaged together or separately

| MODEL | SPECIFICATION (Guage) | NEEDLE LENGTH (mm) | EXTENTION TUBE | COLOR CODE | |
|---------------------|--------------------------|-----------------------|----------------|------------|--|
| 3S-A.D.Cath-16G-45 | | 45 | X | | |
| 3S-A.D.Cath-16G-62 | 16 | 62 | | CDEV | |
| 3S-A.D.Cath-16G-E45 | | 45 | 0 | GREY | |
| 3S-A.D.Cath-16G-E62 | | 62 | U | | |
| 3S-A.D.Cath-18G-45 | | 45 | × | - GREEN | |
| 3S-A.D.Cath-18G-62 | 18 | 62 | ^ | | |
| 3S-A.D.Cath-18G-E45 | 18 - | 45 | | | |
| 3S-A.D.Cath-18G-E62 | | 62 | 0 | | |



- Sharp Needle tip to minimize patient trauma.
- Echogenic needle tip enhances the visualization under ultrasound.
- The design of the handle facilitates the insertion and rotation of the needle during procedure.

| MODEL | SPECIFICATION (Guage) | NEEDLE LENGTH (mm) | TUBE LENGTH | COLOR CODE | |
|-------|--------------------------|--------------------|-------------|------------|--|
| 17G-5 | 17 | 350 | 600 | RED | |
| 17G-6 | 17 | 350 | 900 | NED | |
| 18G-1 | 18 | 350 | 600 | PINK | |
| 18G-2 | 10 | 350 | 900 | T HVIX | |
| 19G-1 | . 19 | 350 | 600 | BROWN | |
| 19G-2 | 19 | 350 | 900 | - BROWN | |

3S-E.T.Cath



- The surface of the tube is very smooth and the tip is manufactured elaborately to prevent injury while injecting to the womb.
- The tube is kink resistance while inject or you may bend with your intended angle.
- Marking for using the adequate length when inserting.

| MODEL | SPECIFICATION (Fr) | CATHETER LENGTH (mm) | COLOR CODE | REMARKS |
|-------------------|-----------------------|----------------------|------------|----------------|
| 3S-E.T.Cath-3.5 | 3.5 | 120 | VIOLET | W/O Depth mark |
| 3S-E.T.Cath-3.5-A | 3.3 | 120 | VIOLET | Depth mark |
| 3S-E.T.Cath-5.5 | 5.5 | 100 | GREEN | W/O Depth mark |

3S-I.U.I.Cath for intrauterine insemination

- Minimize injury while injecting.
- Smooth side hole on the tube.
- Smooth round tip of the tube.
- Smooth surface of the tube.
- Advanced accuracy of operation.
- Easy to handle with intended angle.
- Kink resistance while injecting.

| MODEL | SPECIFICATION (Fr) | CATHETER LENGTH (mm) | COLOR CODE | REMARKS |
|----------------------|-----------------------|----------------------|------------|------------|
| 3S-A.F.Cath-5.5-A100 | 5.5 | 100 | ORANGE | Depth mark |

3S-Endosampler




Biopsy, Used for Biopsy uterine contents.

| MODEL | SPECIFICATION (mm) | CATHETER LENGTH (mm) | COLOR CODE | REMARKS |
|--------------|--------------------|----------------------|------------|-----------------------------|
| Biopsy-EH | | | | End Hole |
| Biopsy-S1-3 | 3 | 245 | WHITE | Side Hole(1ea) |
| Biopsy-S2-3 | . 3 | | VVIIIE | Side Hole(2ea) |
| Biopsy-ES2-3 | | | | End Hole, Side Hole(2ea) |

3S-F.E.Cath

for urinary incontinence



Flexible Endoscopic catheter used to inject or insert drugs into the bladder and urethra using a cystoscope or endoscope.

| MODEL | SPECIFICATION (Fr) | CATHETER LENGTH (mm) | NEEDLE GUAGE (G) | NEEDLE LENGTH (mm) |
|------------------|-----------------------|----------------------|---------------------|--------------------|
| 3S-F.E.Cath-718 | 7 | 400 | 18 | 10 |
| 3S-F.E.Cath-520 | 5 | 350 | 20 | 7 |
| 3S-F.E.Cath-5204 | 5 | 350 | 20 | 4 |

3S-Trocar Set

for spinal access



Catheter, General Purpose.

3S-Trocar Set from Dukwoo Medical is a set of single use sterile radiopaque devices (metal, catheter and needle/trocar) intended to be used for image-guided (using live x-ray), single-puncture, percutaneous access to facilitate an interventional radiology disc decompression procedure.

- Part used: Backbone, cervical vertebral or shoulder, leg joint of the patients.

Needle Tip: Tri facet Trocar or Blunt Tip

| MODEL | SPECIFICATION (Fr / G) | CATHETER LENGTH | REMARK (Outer Tube) |
|----------------------|---------------------------|-----------------|------------------------|
| PDN1010 Cath-6Fr-80 | 6Fr. | 80 | ETFE |
| PDN1010 Cath-6Fr-150 | G11. | 150 | |
| PDN1020 Cath-16G-80 | 16G. | 80 | SUS |
| PDN1020 Cath-16G-150 | 100. | 150 | 303 |
| 3S-16G-100 | 166 | 100 | ETFE |
| 3S-16G-045 | 16G. | 45 | |
| PDN3020 Cath-16G-80 | 166 | 80 | SUS |
| PDN3020 Cath-16G-150 | 16G. | 150 | 303 |
| DM4020 Cath-14G-80 | 446 | 80 | CHC |
| DM4020 Cath-14G-150 | 14G. | 150 | SUS |



| MODEL | SPECIFICATION (Fr / G) | CATHETER LENGTH (mm) | REMARKS (Outer Tube) |
|----------------------|---------------------------|----------------------|-------------------------|
| 14G-80 | 14G. | 80 | PP |
| 3S-G.P.Cath-16G-S300 | 16G. | 300 | SUS |
| 3S-G.P.Cath-5.5-A100 | 5.5Fr. | 100 | PP |
| 3S-G.P.Cath-5-80 | 5Fr. | 80 | PP |
| 3S-G.P.Cath-5Fr-80 | 5Fr. | 80 | PP |
| 5-75 | 5Fr. | 75 | PP |
| 5Fr-75 | 5Fr. | 75 | PP |
| 127G-S25 | 27G | 25 | SUS |

Facilities We have acquired ISO which is an international certification, and is sparing no efforts to establish the optimized quality management system by continuously improving the organization, measuring performance and measuring product quality.




Clean room System



Automation System



Tip Forming



Band Swaging



EO sterilizer



Auto Vision Inspection



Dukwoo Medical Co., Ltd.

HOMEPAGE

http://www.dukwooco.co.kr

TEL

+82-31-355-9330

FAX

+82-31-355-9335

E-MAIL

idukwoo@daum.net

ADDRESS

35, Hwaseong-ro 1616Beongil, Bibong-myeon, Hwaseong-si, Gyeonggi-do, Korea 18284